

Size Information

Size	Women Shoe	Men Shoe
XS	3 1/2 - 5 1/2	2 - 4
SML	6 - 8	4 1/2 - 7
MED	8 1/2 - 11 1/2	7 1/2 - 10 1/2
LRG	11 1/2 - 13 1/2	10 1/2 - 12 1/2
XL	13 1/2 +	12 1/2 +

All walkers are available for 11" and 17"

Washing Instruction

Wash with cold water, gentle cycle, mild soap.
Do not dry clean or tumble dry.
Do not use bleach.
Dry in an air current.
Do not expose to heat source.

Caution

Do not water walker if warmth or heat causes allergic reactions. consult your physician.



Application Instruction

