

14mm 16mm 18mm 20mm 22mm 24mm 26mm 28mm 30mm




HEALTHCARE **HOME MEDICAL SUPPLY**
24 HR IMMEDIATE RESPONSE

Flange Sizing Guide

Print this sheet at 100% custom scale for accuracy. Cut along the blue dotted lines to measure your nipple diameter.

Where to Measure

Your flange size should be 3-5mm larger than your nipple diameter.



*Be sure to measure BEFORE pumping.

Helpful Hint

If you do not have access to a printer, use a coin to determine your nipple diameter.



32mm

34mm

36mm

If your nipple size is...

your flange should measure...

15-17mm

20mm

19-21mm


24mm

23-25mm

28mm

27-29mm

32mm



Call or Visit Our Website

Toll Free
877-240-7 DME

www.healthcaredme.com

For Further Details Please Contact Us.
 Address: 2911 Carpenter Rd, Ann Arbor, MI 48108.
 Phone: (734) 975-6668, Fax: (734) 975-6678
 Email: info@healthcaredme.com

DmeHealthcare HealthcareHomeMedical