


2920 Centennial Rd., Toledo OH 43617-1833 toll-free:

Purchase Order No Date			Fitter's Phone No.	Fitter's Phone No.			
Ordered by			Patient		Age		
Street			Ship to				
City	State	Zip	Street				
Fitter			City	State	Zip		

### **Letters** (Sleeves & Gauntlets)

- A. 3" from wrist (O)
- B. 6" from wrist (O)
- C. Elbow
- Halfway between C and E
- E. Axilla
- Palm at web of thumb
- Palm at base of fingers
- Between thumb joints
- O. Wrist, ahead of wrist bone toward hand



**Optional Shoulder Strap Measurements** (1, 2 and 3)

_EI	Э с	ВА	0 / I	
		4		
		6"	F	G

## **SLEEVES**

LEFT				RIGHT				
	Girth	Length			Girth L		ength	
Ε		O to E		Е		O to E		
D		O to D		D		O to D		
С		O to C		С		O to C		
В		O to B	6" (15.2 cm)	В		O to B	6" (15.2 cm)	
Α		O to A	3" (7.6 cm)	Α		O to A	3" (7.6 cm)	
0				0				

office:


+1 419 474-2973

+1 800 537-1968

## GAUNTLETS ☐ Separate ☐ One Piece with Sleeve

LEFT				RIGHT			
	Girth	Length	l		Girth	Lengt	th
F		O to F		F		O to F	
G		O to G		G		O to G	
Н				Н			
0				0			

# Form F-036 IMF - for Use Without Paper Tapes



Please refer to our current price list for a description of these items

#### **□ HOTLINE SERVICE**

(Guaranteed delivery in 4-5 business days, or less) Additional Cost

## **REQUIRED Allergies:**

## **OPTIONAL SHOULDER STRAP MEASURMENTS:** (See diagram)

- 1. Length from end of last measurement to base of neck, with arm at side (same side as sleeve or gauntlet)
- 2. Waist Circumference
- 3. Diagonal Length for Shoulder Flap with retaining strap

## **OPTIONS & EXTRAS:**

- Zipper: □ Lateral (Outside) (or) ☐ Medial (Inside) (or) ☐ Posterior (Top)
- Garment Color No:
- Elbow Lining: ☐ Inside (or) ☐ Full (or) ☐ Outside
- Strap:
  - Worn to opposite axilla (or)
  - ☐ Adjustable, axilla to waist (or)
  - ☐ Two piece, with separate belt
- Flex/Contracture Seam: □ Elbow
  - Elastic Band: Regular □ 1" □ 2" (or) Silicone □ 1" □ 2"

**DIAGNOSIS:** (Must be indicated to obtain proper counter-pressure)

- □ 22-28 mm Hg Burns
- □ 25-30 mm Hg Lymphedema Mild
- □ 30-35 mm Hg Lymphedema Moderate
- □ 35-40 mm Hg Lymphedema Severe
- ☐ Check this box if there are instructions on other side